

Evolve[®]

Carbomer 980

EASY SQUEEZE BOTTLE

Soft bottle to improve ease of use



23379104

Carbomer 980 0.2% w/w Eye Gel

Preservative-free

10g

Instructions for use

It is important to read this leaflet carefully before using Evolve[®] Carbomer 980. The information provided summarises the important points about the product.



1. What is Evolve[®] Carbomer 980?

Evolve[®] Carbomer 980 is an eye gel containing Carbomer 980 0.2% w/w in a 10g preservative-free multi-dose bottle.

2. Composition

Evolve[®] Carbomer 980 contains Carbomer 980 0.2% w/w, Sodium Chloride, Sorbitol, pH modifier (Sodium Hydroxide) and Water for injections.

3. What is Evolve[®] Carbomer 980 used for?

Evolve[®] Carbomer 980 is a lubricating eye gel intended for use in the treatment of Dry Eye Disease (Dry Eye Syndrome), where it helps by providing additional lubrication and protection to the surface of your eyes. Evolve[®] Carbomer 980 is **preservative-free, phosphate-free, buffer-free and ideal for sensitive eyes**. The Evolve[®] system has been designed to combine the benefits of unit dose preservative-free eye gel with the convenience of a standard eye drop bottle.

Evolve[®] Carbomer 980 is compatible with contact lenses.

The bottle delivers 360 doses of solution.

For use in adults. Clinical data to support use in persons <18 years is not available.

The bottle has been designed to assist ease of use for patients. The body of the bottle is easy to grip using a pinch grip between the thumb and forefinger and has been constructed of a soft pliable plastic making squeezing and control of dispensing a gel easy and convenient. Only GENTLE pressure is required to squeeze one gel drop out. The blue tip contrasts against a white background, this will assist in aiming the gel drop as you invert the bottle towards your eye. Lastly, there is no protruding spout that could poke and damage your eye.

4. How to use

Before opening a new bottle of Evolve[®] Carbomer 980, please check the safety seal of the bottle is unbroken by checking the integrity of the six bridges between the tamper-evident ring and the cap.

- Wash and dry your hands well before use.
- Rotate the cap anticlockwise to break the seal.
- Remove the cap.
- Tilt your head back and look up.
- Hold the bottle between the thumb and forefinger of your hand.
- Gently pull down the lower lid of the eye to form a small pocket between your eyelid and your eye.
- Apply moderate pressure in the middle of the bottle to release one drop of Evolve[®] Carbomer 980 into the pocket. Please note that there may be a few seconds delay between squeezing and the drop coming out. **Do not squeeze too hard.**
- Blink a few times to spread the drop over the eye.
- If you think that you have missed the eye, then insert another drop.
- Carefully wipe away any excess product around the eyelid with a tissue.
- Repeat the process for the other eye, if required.
- Shake the bottle once downwards to remove any excess drop (see Figure 1) before replacing the cap on the container, do not touch the applicator tip with anything, including your eye or your fingers.



Figure 1

If you are unsure about how to use **Evolve® Carbomer 980**, consult your healthcare professional.

5. How much to use?

The usual dosage is one drop to be applied two times a day in each eye.

6. Precautions and Warnings

- For topical ocular use only. Do not use if the safety seal on the cap of the bottle is broken.
- Do not use if you are allergic to any of the ingredients.
- Do not touch the tip of the bottle to your eye or any other surface.
- Do not share this product with any other person.
- Keep out of the reach and sight of children.
- Do not use after the expiry date shown on the carton or printed on the bottle.
- This product can cause transient blurring of vision after application. Wait until your vision is clear before driving or using machinery.
- Leave at least 5 minutes between administrations of any other eye drops or eye ointments. For use of **Evolve® Carbomer 980** with contact lenses, please consult your healthcare practitioner.
- Do not use if you are pregnant, breast-feeding or have recently undergone an ophthalmic procedure.
- Medicom Healthcare Ltd. does not advise use of this product in a way that is different from its label.
- Your healthcare professional may advise you to use this product in a different way, so follow their directions and read the label.
- If your eye condition persists, gets worse or you experience any unusual side effects after using this product, please consult your healthcare professional and report to Medicom Healthcare Ltd.

7. Contraindications

Do not use if you are allergic to Carbomer or to any of the ingredients listed in section 2-Composition.

8. Possible Side Effects

- Transient blurring of vision
- Mild discomfort in eye
- Possible allergic reaction

If you experience any additional adverse events, consult your healthcare professional and report to Medicom Healthcare Ltd.

9. How to store Evolve® Carbomer 980?

Do not store above 25°C.

Discard the **Evolve® Carbomer 980** six months after first opening.

Dispose of the bottle according to your local waste management regulations.

Bottle Discards Date Reminder: Patients using this product are advised to fill out their name, first open date & discard date in the designated area on the top inside flap of the carton box. This will enable the patient to know exactly when to discard the product after first opening.

This leaflet applies to **Evolve® Carbomer 980** and not to any other product.

Last revision of the text: Aug 2025



Discard the contents 6 months after first opening



Do not use if package is damaged and consult instructions for use



Sterilised using aseptic processing techniques



Consult instructions for use



Do not store above 25°C. Do not freeze the product



Medical device symbol



Medicom Healthcare Ltd.,
Lynton House,
7-12 Tavistock Square, Kings Cross,
London, WC1H 9LT, U.K.

**UK
CA
0120**