

20 Individually Sealed Sachets. Sterile and Detergent Free Wipes.



- ✓ For adults & babies
- ✓ Gently Cleans away debris caused by ocular infection and allergy
- ✓ Daily eyelid hygiene
- ✓ Cooling & Soothing

For MGD, blepharitis and all types of dry eye

Evolve® Eyelid Wipes

x20

10-ELWUKCA-21-GB-CV1

Contains: Sodium Chloride, Potassium Chloride, Calcium Chloride Dihydrate, Magnesium Chloride Hexahydrate, Sodium Acetate, Sodium Citrate Dihydrate, Water Purified.

Directions

Before using a new wipe, please check the integrity of the sachet. Do not use if it is opened or damaged.

1. Wash and dry your hands well before use.
2. Remove wipe from sachet.
3. If you are a contact lens wearer remove your lenses before use.
4. Close the eye and wipe from nasal corner of eye outwards across the lashes.
5. After use, wipes to be disposed of immediately.
6. Use a new wipe for each eye



PIP CODE: 404-7577



UK
CA
0120

MD

Medicom Healthcare Ltd.,
Lynton House,
7-12 Tavistock Square,
Kings Cross, London
WC1H 9LT, U.K.

Evolve® Eyelid Wipes provide soothing relief for the symptoms of dry eye disease, including tired, sore, gritty, and watery eyes.

The symptoms of dry eye are often caused by problems with the glands in your eyelids (Meibomian Gland Dysfunction or MGD) or infection in your eyelids (blepharitis).

Evolve® Eyelid Wipes have been specifically designed to tackle these problems and are suitable for adults and babies. The textured wipe is designed to remove any debris around the eyelids and eyelashes helping restore good eyelid hygiene.

Evolve® Eyelid Wipes are preservative-free and contain no detergents. Detergents can destabilise your tears which can lead to dry eye. **Evolve® Eyelid Wipes** contain water and electrolytes which mimic your natural tears

Evolve® Eyelid Wipes are ideal for daily use to maintain eyelid hygiene as part of the treatment of blepharitis and for cleaning discharge or crustiness from the eyelid and eyelashes because of ocular allergy or ocular infections. These types of products are recommended by eye care professionals.

Evolve® Eyelid Wipes

x20

For MGD, blepharitis and all types of dry eye

- ✓ Cooling & Soothing
- ✓ Daily eyelid hygiene
- ✓ Gently Cleans away debris caused by ocular infection and allergy
- ✓ For adults & babies



20 Individually Sealed Sachets.
Sterile and Detergent Free Wipes.

Precautions and Warnings

- Do not use if the sachet is open or damaged.
- Do not use if you are allergic to any of the ingredients.
- Do not reuse.
- Single patient use only.
- Sterile until open.
- Avoid contact with the eye surface.
- For external use only, do not swallow.
- Do not use after the expiry date shown on the carton and punched on the sachet.
- Keep out of reach and sight of the children.
- Dispose of the used wipe and sachet according to your local waste management regulations.
- Store between 1°C and 25°C.
- If your eye condition persists consult your healthcare professional.

NVZ